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Instruction for use Tener®IT

Intended use

Tener®IT is a spring hard, beryllium- and nickel-free, cobalt based alloy for the manufacture of partial dentures.

General guidelines for handling

This instruction for use includes important processing steps and recommendations for **Tener®IT**.

Indication

For the production of model cast works for removable restorations.

Design

The design has to be made on an investment model (e.g. **Tener®Vest**, Eisenbacher Dentalwaren ED GmbH) under consideration of the standard rules of designs for partial dentures. The modelation should be done with wax that fire without leaving residues.

Sprue design

The sprues should have a diameter of 3.5 mm with a length of 20-30 mm. Massive areas should additionally have a melt reservoir (ø 5mm). The cast hopper should be 4mm above the highest point of the partial denture wax design.

Investing

Mixing ratios, expansion control, preheating processes are different depending on the investment material. Please refer to the processing guide of the corresponding investment material (e.g. **Tener®Vest**, Eisenbacher Dentalwaren ED GmbH). Recommended pre-heating temperature: 900 – 950°C.

Melting and casting

Tener®IT should be melted in a ceramic crucible. Please do not use graphite crucibles and no flux! Avoid the overheating of the melt. Prevent multiple casts of melt bottoms. The chemical and mechanical properties can only be guaranteed for new material. Melting with open flame (acetylene / oxygen) and inductive melting: Once the cylinders are melted and a so-called shadow falls across the molten metal, before the oxide skin begins to split, start the casting. **Maximum temperature for casting: 1504 °C.**

Devesting

Let the muffle cool down to room temperature (ca. 20°C), do not quench with water. Put the cooled muffle into water to avoid dust generation during the devesting. Sandblast the surface with 250 µm of aluminium oxide with 3-4 bar. Clean the partial denture base with a steam cleaner.

Elaboration

The framework can be elaborated with standard carbide cutters or aluminium oxide stones and rubbers, look for smooth transitions.

Soldering / Laser welding

For soldering we recommend **Kera®Wel-con**, Eisenbacher Dentalwaren ED GmbH. The **Tener®IT** alloy is also approved for laser-welding, we recommend **Kera®LA-con**, Eisenbacher Dentalwaren ED GmbH.

Handling conditions / Safety

Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.

Contraindications and side effects

If the instructions are observed during the production processes, incompatibilities with CoCr alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the contra-indications and side effects.

Disposal Instructions

Consult the material safety data sheets or national regulations for disposal.

Storage conditions

Temperature, humidity or light has no effect on the product properties.

Quantity

Please consider the label on the package.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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